

Invitation to EWMA conference

We would like to invite you to the European Wound Management Association Conference, which will be held from the 13th to 15th May in London. At the expert event for healthcare professionals focussing on various aspects of treatment for wounds and pressure ulcers, the LINET Group company will present a wide range of innovative solutions aimed at the prevention and treatment of pressure ulcers.



It will be our great pleasure to welcome you at stand 2-70, where you can find out about various technologies intended for intensive and long-term care. **We look forward to meeting you!**

Multicare + Symbioso 200, Eleganza3XC + Virtuoso + Sentida 7i + Air 2Care

Two professionals: Multicare and Symbioso

The **Symbioso mattress integrated in the Multicare resuscitation bed** is exceptionally silent and effective. The mattress works on the **dual principle of maintaining a constant pressure and suitable microclimate**. The combination of both systems provides effective prevention against the formation of pressure ulcers by eliminating pressure and dampness between the patient's body and the mattress. A suitable microclimatic effect is ensured by the constant airflow in the mattress. In combination with the **positioning options of the Multicare bed**, it has a significant preventative and therapeutic effect.

Virtuoso

You can see the Eleganza 3XC intensive care bed at stand 2-70 in combination with the Virtuoso active anti-decubitus system. This mattress works on the principle of **alternating 3 linked cells**, which provide **relief and reduction in pressure to zero on the patient's entire body every 7.5 minutes**. The Virtuoso mattress delivers an effective **24-hour care** solution for patients with the greatest risk of pressure ulcer formation.

Intelligent solution for long-term care

The **new Sentida 7i bed** builds on the previous models from many aspects, while maintaining their practical and valued functions - **the low height of the mattress platform and multifunctional side rails** which, in addition to their safety function, support the **mobilisation of clients**. The intelligent technologies are innovative – for example **control panel in the form of a touchscreen display** and safety sensors, such as the **wet bed sensor**. Dampness is a negative factor influencing the formation of pressure ulcers. The personnel can react immediately and **eliminate potential danger and discomfort** thanks to timely detection. Besides the bed, the new products include the Air2Care mattress, also intended for **patients who remain recumbent** for a long time. The **variability of use** is one of the greatest advantages of this mattress. While it can be used on **top of an existing mattress**, there is also a **self-supporting system** available with a bottom foam layer or a static air layer. **The air layer** of the active mattress **consists of a structure of joined cells**, which are inflated alternately **at 12-minute intervals**. The mattress thus ensures the **regular reduction of pressure** as one of the most dangerous factors leading to the formation of pressure ulcers.



EWMA 2015 (EUROPEAN WOUND MANAGEMENT ASSOCIATION)

13-15 May, 2015

Excel, London

STAND NUMBER 2-70